

# Costech 4010 Elemental Analyzer Sample Size Determination v0.2

Robin B. Trayler

---

Please contact the Stable Isotope Ecology Laboratory Technical Director with any questions or clarifications. **StableIsotopes@ucmerced.edu**

---

## Introduction

Determining the correct sample size for  $\delta^{13}\text{C}$  and  $\delta^{15}\text{N}$  analyses is important because the mass spectrometer signal intensity, sensitivity, and reproducibility is related to sample mass. These notes make three assumptions. 1. The the carbon and nitrogen content of a material is known prior to analysis. 2. Analysis is taking place in the University of California, Merced Stable Isotope Ecology laboratory using a *Costech 4010 Elemental Analyzer*.

\* Other laboratories have similar instrumentation setups, however instrument sensitivity will vary between configurations.

3. The *Delta V+ IRMS*  $\text{CO}_2$  sample dilution is set to 89%. \* This is the default configuration for routine analyses. Other configurations are possible but you should contact the Technical Director prior to sample preparation.

## Background

The signal intensity of a sample gas measured by a *Delta V+ isotope ratio mass spectrometer* is directly related to the amount of gas. In the case of the *Costech 4010 Elemental Analyzer* the signal intensity is linearly related to sample mass (Figure 1).




Figure 1: Sample carbon and nitrogen content versus signal intensity.

The *Delta V+ isotope ratio mass spectrometer* has a measurement range from 0 to 50

24 volts. However the sensitivity is best below about 30 volts. Furthermore, both  $\delta^{13}\text{C}$  and  
 25  $\delta^{15}\text{N}$  have have a non-linear relationship to signal intensity. This means that a sample  
 26 with a small signal intensity will have a different measured  $\delta$  value than a sample of the  
 27 same material at a large signal intensity.




Figure 2: Mass dependence of  $\delta^{13}\text{C}$  and  $\delta^{15}\text{N}$  values.

28 The  $\delta$ -value-mass-dependence is linear and lowest between about 5 to 15 volts for  $\delta^{13}\text{C}$   
 29 and 5 to 25 for  $\delta^{15}\text{N}$  (Figure 2). Since the relationship between signal intensity (volts) is

linearly related to carbon and nitrogen content (Figure 1), this equates to the mass range shown in Table 1.

Table 1: minimum and maximum carbon and nitrogen contents for  $\delta^{13}\text{C}$  and  $\delta^{15}\text{N}$  analysis

isotope	minimum (mg)	maximum (mg)
$\delta^{13}\text{C}$	0.2	0.8
$\delta^{15}\text{N}$	0.05	0.2

## Calculating Sample Size

If the carbon and nitrogen content of a material is known prior to analysis you can calculate the optimal sample size using the following equations to calculate the range of sample weights for carbon and nitrogen.

### $\delta^{13}\text{C}$ Range

$$Sample\ size_{min} = 0.2 \div (sample\ carbon\% \div 100)$$

$$Sample\ size_{max} = 0.8 \div (sample\ carbon\% \div 100)$$

### $\delta^{15}\text{N}$ Range

$$Sample\ size_{min} = 0.05 \div (sample\ nitrogen\% \div 100)$$

$$Sample\ size_{max} = 0.2 \div (sample\ nitrogen\% \div 100)$$

## Measuring both $\delta^{13}\text{C}$ and $\delta^{15}\text{N}$

It is difficult to concurrently measure  $\delta^{13}\text{C}$  and  $\delta^{15}\text{N}$  in samples with high C:N ratios. Based on the sample ranges given in Table 1 the maximum C:N ratio is about 16 (*carbon maximum*  $\div$  *nitrogen minimum* or  $0.8 \div 0.05$ ) and the minimum C:N ratio is 1 (*carbon minimum*  $\div$  *nitrogen maximum* or  $0.2 \div 0.2$ ). For analyses of materials with C:N ratios outside this range please consult the Technical Director prior to sample preparation.

If your samples have a C:N ratio between about 1 to 16 then  $\delta^{13}\text{C}$  and  $\delta^{15}\text{N}$  can be measured concurrently. Calculate the range of acceptable sample sizes using the equations for carbon and nitrogen shown above and then choose a final sample size that overlaps both results. If your samples have a C:N ratio significantly outside this range then  $\delta^{13}\text{C}$  and  $\delta^{15}\text{N}$  may need to be measured separately. Please consult the Technical Director for guidance.

## Commonly Analyzed Materials

Table 2: Carbon and nitrogen contents of commonly analyzed materials.

material	carbon %	nitrogen %
acetanilide	71.09	10.36
biochar	> 80	< 1.5
collagen	40 - 45	10 - 15
glutamic acid	40.3	9.52
keratin	40 - 45	10 - 15
plants	40 - 50	< 1 - 5

material	carbon %	nitrogen %
proteins	40 - 45	10 - 15
soil	highly variable	highly variable

## Worked Example

Acetanilide is common reference material that is 71.09% carbon and 10.36% nitrogen (Table 2).

### $\delta^{13}\text{C}$

$$sample\ size_{min} = 0.28\ mg = 0.2 \div (71.09\% \div 100)$$

$$sample\ size_{max} = 1.13\ mg = 0.8 \div (71.09\% \div 100)$$

### $\delta^{15}\text{N}$

$$Sample\ size_{min} = 0.48\ mg = 0.05 \div (10.36\% \div 100)$$

$$Sample\ size_{max} = 1.93\ mg = 0.2 \div (10.36\% \div 100)$$

In this case the minimum sample size for  $\delta^{13}\text{C}$  measurements is lower than the minimum for measuring  $\delta^{15}\text{N}$ . the minimum sample size for  $\delta^{15}\text{N}$  analysis is in the middle of the range for  $\delta^{13}\text{C}$  so a sample size of **0.5 mg** would be appropriate.